

ARG10117 anti-Fluorescein antibody [2A3]

Package: 500 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [2A3] recognizes Fluorescein
Tested Reactivity	Other
Tested Application	IHC-Fr, IHC-P
Specificity	This antibody reacts with both free and conjugated FITC.
Host	Mouse
Clonality	Monoclonal
Clone	2A3
Isotype	IgG2a
Target Name	Fluorescein
Species	Others
Immunogen	FITC-BSA conjugate
Conjugation	Un-conjugated

Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--------------------------------------------------------------------------------------------------------------------------------------------

Properties

Form	Liquid
Purification	Protein A affinity purified.
Buffer	PBS (pH 7.4) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Concentration	1.0-2.0 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Controls and Markers antibody
---------------	-------------------------------