

## ARG10163 anti-Triiodothyronine antibody [3A6]

Package: 100 µg, 50 µg  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [3A6] recognizes Triiodothyronine
Tested Reactivity	Other
Tested Application	ELISA
Host	Mouse
Clonality	Monoclonal
Clone	3A6
Isotype	IgG1
Target Name	Triiodothyronine
Antigen Species	Others
Immunogen	Conjugate of purified triiodothyronine (Sigma T-2752) with BSA
Conjugation	Un-conjugated

### Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

### Properties

Form	Liquid
Purification	Protein A affinity purified.
Buffer	PBS (pH 7.4) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Concentration	1.0-2.0 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Signaling Transduction antibody
---------------	---------------------------------