

ARG10263 anti-Ricinus communis agglutinin 60 antibody [CP75]

Package: 250 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [CP75] recognizes Ricinus communis agglutinin 60
Tested Reactivity	Plnt
Tested Application	ELISA, WB
Host	Mouse
Clonality	Monoclonal
Clone	CP75
Isotype	IgG1
Target Name	Ricinus communis agglutinin 60
Species	Plant
Immunogen	Full length natural RCA60 protein
Conjugation	Un-conjugated

Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Purification	Protein G affinity purified
Buffer	PBS (pH 7.4) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Concentration	1.0-2.0 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Cell Biology and Cellular Response antibody
---------------	---