

ARG10717 anti-EosFP antibody

Package: 50 µl
Store at: -20°C

Summary


Product Description	Rabbit Polyclonal antibody recognizes EosFP
Tested Reactivity	Other
Tested Application	ICC/IF, IHC-Fr, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	EosFP
Immunogen	Full length recombinant protein based on the Clontech enhanced GFP.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ICC/IF	1:5000
	IHC-Fr	1:5000
	WB	1:1000 - 1:5000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


Properties

Form	Liquid
Purification	Unpurified.
Buffer	Serum and 5mM sodium azide.
Preservative	5mM sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG10717 anti-EosFP antibody ICC/IF image

Immunocytochemistry: Transfected HEK293 cells which overexpress EosFP protein were stained with ARG10717 anti-EosFP antibody. Cells which are transfected with EosFP (left panel) are bright green. On staining with ARG10717 anti-EosFP antibody antibody in red (middle panel), EosFP expressing cells appear orange (right panel). Most HEK293 cells are not transfected so only the nucleus of these cells can be visualized with a blue DNA stain.



ARG10717 anti-EosFP antibody WB image

Western blot: 1) Non-transfected HEK293 crude homogenate, and 2) Crude homogenate of transfected HEK293 cells which overexpress protein EosFP. ARG10717 anti-EosFP antibody at 1: 1000 dilution.
