

## ARG10725 anti-Parvalbumin antibody [3C9]

Package: 50 µl  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [3C9] recognizes Parvalbumin
Tested Reactivity	Hu, Ms, Rat, Cow
Tested Application	ICC/IF, IHC-Fr, WB
Host	Mouse
Clonality	Monoclonal
Clone	3C9
Isotype	IgG1
Target Name	Parvalbumin
Antigen Species	Human
Immunogen	Full-length recombinant Human Parvalbumin.
Conjugation	Un-conjugated
Alternate Names	D22S749; Parvalbumin alpha

### Application Instructions

Application table	Application	Dilution
	ICC/IF	1:1000 - 1:5000
	IHC-Fr	1:1000 - 1:5000
	WB	1:1000 - 1:5000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Calculated Mw	12 kDa	


### Properties

Form	Liquid
Purification	Affinity purification.
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

## Bioinformation


Gene Symbol	PVALB
Gene Full Name	parvalbumin
Background	The protein encoded by this gene is a high affinity calcium ion-binding protein that is structurally and functionally similar to calmodulin and troponin C. The encoded protein is thought to be involved in muscle relaxation. [provided by RefSeq, Jul 2008]
Function	In muscle, parvalbumin is thought to be involved in relaxation after contraction. It binds two calcium ions. [UniProt]

## Images



ARG10725 anti-Parvalbumin antibody [3C9] ICC/IF image


Immunocytochemistry: Adult Rat cerebellum floating section was stained with ARG10725 anti-Parvalbumin antibody [3C9] at 1:1000 (green). Parvalbumin is prominently expressed in the dendrites and perikarya of Purkinje cells and some interneurons in the molecular layer. Blue is a DNA stain.



ARG10725 anti-Parvalbumin antibody [3C9] IHC-Fr image


Immunohistochemistry: Frozen section of Rat cerebellum stained with ARG10725 anti-Parvalbumin antibody [3C9] (green) at 1:1000 dilution and costained with Chicken pAb to calbindin (red) at 1:2000 dilution. DAPI (blue) for nuclear staining. (Sample preparation: Following transcardial perfusion of Rat with 4% paraformaldehyde, brain was post fixed for 24 hours, cut to 45  $\mu$ M, and free-floating sections were stained with above antibodies.)

Most Purkinje cells coexpress both parvalbumin and calbindin and so appear yellow, whereas basket, stellate and Golgi cells express parvalbumin only and so appear green.



ARG10725 anti-Parvalbumin antibody [3C9] WB image

Western blot: Rat skeletal muscle lysate (left) and His-tagged recombinant proteins (right): 1) parvalbumin, 2) calretinin, and 3) calbindin. Blots were stained with ARG10725 anti-Parvalbumin antibody [3C9] at 1:1000. In skeletal muscle lysates, this antibody recognizes a band at 12 kDa which represents parvalbumin. It reacts to recombinant parvalbumin only, not the other calcium-binding proteins listed. Note, this antibody does not recognize parvalbumin in Rat or mouse brain lysates on western blots. Protein from Rat skeletal muscle lysates were transferred to PVDF membrane.



#### ARG10725 anti-Parvalbumin antibody [3C9] IHC-Fr image

Immunohistochemistry: Frozen section of Rat brain cerebellum stained with ARG10725 anti-Parvalbumin antibody [3C9] (green) at 1:1000 dilution and costained with Chicken pAb to calbindin (red) at 1:2000 dilution. (Sample preparation: Following transcardial perfusion of Rat with 4% paraformaldehyde, brain was post fixed for 24 hours, cut to 45  $\mu$ M, and free-floating sections were stained with above antibodies.)

Most Purkinje cells coexpress both parvalbumin and calbindin and so appear yellow, whereas basket, stellate and Golgi cells express parvalbumin only and so appear green.