

ARG10742 anti-Nop2p antibody [22G1]

Package: 100 µl
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [22G1] recognizes Nop2p
Tested Reactivity	Yeast
Tested Application	ICC/IF, IHC-Fr, WB
Host	Mouse
Clonality	Monoclonal
Clone	22G1
Isotype	IgG3
Target Name	Nop2p
Antigen Species	Yeast
Immunogen	Yeast nuclear preparations.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ICC/IF	1:1000 - 1:5000
	IHC-Fr	1:1000 - 1:5000
	WB	1:2000 (cell lysates), 1:10000 (nuclear fractions)
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Unpurified.
Buffer	Tissue culture supernatant.
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.