

**ARG10745**  
**anti-Nsr1p antibody [31C4]**Package: 100 µl  
Store at: -20°C

### Summary


Product Description	Mouse Monoclonal antibody [31C4] recognizes Nsr1p
Tested Reactivity	Yeast
Tested Application	ICC/IF, WB
Host	Mouse
Clonality	Monoclonal
Clone	31C4
Isotype	IgG
Target Name	Nsr1p
Species	Yeast
Immunogen	Yeast nuclear preparations.
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ICC/IF	1:500 - 1:5000
	WB	1:1000 - 1:10000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG10745 anti-Nsr1p antibody [31C4] WB image

Western blot: Whole yeast protein extracts with a collection of arigo's antibodies. The blot for ARG10745 anti-Nsr1p antibody [31C4] is in the indicated lane (clone name), and the number indicates the SDS-PAGE molecular weight in kDa.