

ARG11078 anti-Brucella abortus antibody [BrF11]

Package: 200 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [BrF11] recognizes Brucella abortus
Tested Reactivity	Bacteria
Tested Application	ELISA, WB
Specificity	LPS of Brucella abortus and with intact Brucella-cells.
Host	Mouse
Clonality	Monoclonal
Clone	BrF11
Isotype	IgG1
Target Name	Brucella abortus
Species	Bacteria
Immunogen	Brucella abortus
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purification with Protein G.
Buffer	PBS (pH 7.4) and 0.09% Sodium azide.
Preservative	0.09% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.