

## ARG11090 anti-MERS-CoV nucleocapsid protein antibody [7H6]

Package: 100 µg  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [7H6] recognizes MERS-CoV nucleocapsid protein
Tested Reactivity	Virus
Tested Application	ELISA, ICC/IF, IP, WB
Host	Mouse
Clonality	Monoclonal
Clone	7H6
Isotype	IgG2b
Target Name	MERS-CoV nucleocapsid protein
Species	Virus
Immunogen	GST-tagged MERS-CoV nucleocapsid protein corresponding to aa. 196-413.
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	ICC/IF	Assay-dependent
	IP	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Purification with Protein A.
Buffer	PBS and 0.02% Sodium azide.
Preservative	0.02% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

## Bioinformation

---

**Gene Full Name** Middle East Respiratory Syndrome coronavirus nucleocapsid protein

**Highlight** Related products:  
[MERS antibodies](#); [Anti-Mouse IgG secondary antibodies](#);  
Related news:  
[ACE2, receptor of 2019-nCoV](#)