

ARG20592 Mouse IgM Isotype Control antibody [11E10] (low endotoxin)

Package: 100 µg
Store at: -20°C

Summary

Product Description	Azide free and low endotoxin Mouse Monoclonal antibody [11E10] as a negative control antibody for Mouse IgM
Tested Application	BL, ELISA, FACS, FLISA, IHC-Fr, IHC-P
Specificity	Lipopolysaccharide
Host	Mouse
Clonality	Monoclonal
Clone	11E10
Isotype	IgM, lambda
Target Name	Mouse IgM
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	ELISA	< 1 µg/ml
	FACS	< 1 µg/10 ⁶ cells
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Purified
Purification Note	Low endotoxin
Buffer	PBS
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.