

ARG20603 Mouse IgG1 Isotype Control antibody [15H6] (APC-Cyanine 7)

Package: 50 tests
Store at: 4°C

Summary

Product Description	APC-Cyanine 7-conjugated Mouse Monoclonal antibody [15H6] as a negative control antibody for Mouse IgG1
Tested Application	BL, ELISA, FACS, FLISA, ICC/IF, IHC-Fr, IHC-P, WB
Specificity	T-2 mycotoxin
Host	Mouse
Clonality	Monoclonal
Clone	15H6
Isotype	IgG1, kappa
Target Name	Mouse IgG1
Conjugation	APC-Cyanine 7

Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	ELISA	Assay-dependent
	FACS	10 µl/10 ⁶ cells
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.