

ARG20616 Mouse IgG2b Isotype Control antibody [A-1] (low endotoxin)

Package: 100 µg
Store at: -20°C

Summary

Product Description	Azide free and low endotoxin Mouse Monoclonal antibody [A-1] as a negative control antibody for Mouse IgG2b
Tested Application	ELISA, FACS, FLISA, IHC-Fr
Specificity	Chicken IgA
Host	Mouse
Clonality	Monoclonal
Clone	A-1
Isotype	IgG2b, kappa
Target Name	Mouse IgG2b
Target Ig	Chicken IgA
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	1 - 5 µg/ml
	FACS	< 1 µg/10 ⁶ cells
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified
Purification Note	Low endotoxin
Buffer	PBS (pH 7.4), low endotoxin
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.