

ARG20631 Mouse IgA Isotype Control antibody [S107] (FITC)

Package: 50 tests
Store at: 4°C

Summary


Product Description	FITC-conjugated Mouse Monoclonal antibody [S107] as a negative control antibody for Mouse IgA
Tested Application	FACS, FLISA
Specificity	C-carbohydrate
Host	Mouse
Clonality	Monoclonal
Clone	S107
Isotype	IgA, kappa
Target Name	Mouse IgA
Conjugation	FITC

Application Instructions

Application table	Application	Dilution
	FACS	10 µl/10 ⁶ cells
	FLISA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG20631 Mouse IgA Isotype Control antibody [S107] (FITC) FLISA image

FLISA: The plate was coated with Goat anti-Mouse IgG antibody, [ARG21539](#) Goat anti-Mouse IgA antibody (pre-adsorbed) and [ARG21517](#) Goat anti-Mouse IgM antibody (pre-adsorbed). Serially diluted ARG20631 Mouse IgA Isotype Control antibody [S107] (FITC) was captured and fluorescence intensity quantified.