

**ARG20634**  
**Mouse IgA Isotype Control antibody [S107] (APC)**Package: 50 tests  
Store at: 4°C

### Summary

Product Description	APC-conjugated Mouse Monoclonal antibody [S107] as a negative control antibody for Mouse IgA
Tested Application	ELISA, FACS, FLISA
Specificity	C-carbohydrate
Host	Mouse
Clonality	Monoclonal
Clone	S107
Isotype	IgA, kappa
Target Name	Mouse IgA
Conjugation	APC

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	10 µl/10 <sup>6</sup> cells
	FLISA	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.