

ARG20641 Mouse IgG Isotype Control antibody (AP)

Package: 500 µl
Store at: 4°C

Summary


Product Description	AP-conjugated Mouse Polyclonal antibody as a negative control antibody for Mouse IgG
Tested Application	BL, ELISA, FACS, FLISA, IHC-P, WB
Host	Mouse
Clonality	Polyclonal
Isotype	IgG
Target Name	Mouse IgG
Conjugation	AP

Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	ELISA	1:500 - 1:5000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent
	Application Note <p>* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.</p>	

Properties

Form	Liquid
Buffer	50 mM Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl ₂ and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG20641 Mouse IgG Isotype Control antibody (AP) ELISA image

ELISA: The plate was coated with Goat anti-Mouse IgG antibody, [ARG21539](#) Goat anti-Mouse IgA antibody (pre-adsorbed), and [ARG21517](#) Goat anti-Mouse IgM antibody (pre-adsorbed). Serially diluted ARG20641 Mouse IgG Isotype Control antibody (AP) was captured and quantified.