

## ARG20653 Rat IgG Isotype Control antibody (APC)

Package: 50 µg  
Store at: 4°C

### Summary


Product Description	APC-conjugated Rat Polyclonal antibody as a negative control antibody for Rat IgG
Tested Application	BL, FACS, FLISA, ICC/IF, IHC-P
Host	Rat
Clonality	Polyclonal
Isotype	IgG
Target Name	Rat IgG
Conjugation	APC

### Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	FACS	< 0.1 µg/10 <sup>6</sup> cells
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	IHC-P	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG20653 Rat IgG Isotype Control antibody (APC) FLISA image

FLISA: The plate was coated with [ARG21957](#) Goat anti-Rat IgG antibody (pre-adsorbed) and [ARG23828](#) Mouse anti-Rat IgM antibody [M2A1]. Serially diluted ARG20653 Rat IgG Isotype Control antibody (APC) was captured and fluorescence intensity quantified.