

ARG20669 Goat F(ab')₂ IgG Isotype Control antibody (Biotin)

Package: 100 µg
Store at: 4°C

Summary

Product Description	Biotin-conjugated Goat Polyclonal antibody as a negative control antibody for Goat F(ab') ₂ IgG
Tested Application	ELISA, FACS, FLISA
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	Goat Fab'2 IgG
Conjugation	Biotin

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	< 1 µg/10 ⁶ cells
	FLISA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.