

ARG20681
Rabbit F(ab')₂ IgG Isotype Control antibody (PE)

Package: 50 tests

Store at: 4°C

Summary

Product Description	PE-conjugated Rabbit Polyclonal antibody as a negative control antibody for Rabbit F(ab') ₂ IgG
Tested Application	FACS, FLISA
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Rabbit Fab'2 IgG
Conjugation	PE


Application Instructions

Application table	Application	Dilution
	FACS	10 µl/10 ⁶ cells
	FLISA	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG20681 Rabbit F(ab')₂ IgG Isotype Control antibody (PE) FLISA image

FLISA: The plate was coated with Goat anti-Rabbit IgG (Fab) antibody and Goat anti-Rabbit IgG (Fc) antibody. Serially diluted ARG20681 Rabbit F(ab')₂ IgG Isotype Control antibody (PE) was captured and fluorescence intensity quantified.