

## ARG20685 Mouse F(ab')<sub>2</sub> IgG1 Isotype Control antibody [15H6] (FITC)

Package: 50 tests  
Store at: 4°C

### Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [15H6] as a negative control antibody for Mouse F(ab') <sub>2</sub> IgG1
Tested Application	ELISA, FACS, FLISA, IHC-P
Specificity	T-2 mycotoxin
Host	Mouse
Clonality	Monoclonal
Clone	15H6
Isotype	IgG1, kappa
Target Name	Fab'2 IgG1
Conjugation	FITC

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	10 µl/10 <sup>6</sup> cells
	FLISA	Assay-dependent
	IHC-P	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.