

## ARG20691 Rat IgG1 Isotype Control antibody [KLH/G1-2-2] (AP)

Package: 500 µl  
Store at: 4°C

### Summary

Product Description	AP-conjugated Rat Monoclonal antibody [KLH/G1-2-2] as a negative control antibody for Rat IgG1
Tested Application	BL, ELISA, FACS, FLISA, ICC/IF
Specificity	KLH
Host	Rat
Clonality	Monoclonal
Clone	KLH/G1-2-2
Isotype	IgG1, kappa
Target Name	Rat IgG1
Conjugation	AP

### Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	ELISA	1:2000 - 1:4000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	50 mM Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl <sub>2</sub> and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.