

ARG20698 Rat IgG1 Isotype Control antibody [KLH/G1-2-2] (PE-Cyanine 7)

Package: 50 µg
Store at: 4°C

Summary

Product Description	PE-Cyanine 7-conjugated Rat Monoclonal antibody [KLH/G1-2-2] as a negative control antibody for Rat IgG1
Tested Application	BL, ELISA, FACS, FLISA, ICC/IF
Specificity	KLH
Host	Rat
Clonality	Monoclonal
Clone	KLH/G1-2-2
Isotype	IgG1, kappa
Target Name	Rat IgG1
Conjugation	PE-Cyanine 7

Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	ELISA	Assay-dependent
	FACS	< 0.2 µg/10 ⁶ cells
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.