

## ARG20702 Rat IgG2a Isotype Control antibody [KLH/G2a-1-1] (FITC)

Package: 100 µg  
Store at: 4°C

### Summary


Product Description	FITC-conjugated Rat Monoclonal antibody [KLH/G2a-1-1] as a negative control antibody for Rat IgG2a
Tested Application	BL, FACS, FLISA, IHC-Fr, IHC-P
Specificity	KLH
Host	Rat
Clonality	Monoclonal
Clone	KLH/G2a-1-1
Isotype	IgG2a, kappa
Target Name	Rat IgG2a
Conjugation	FITC

### Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	FACS	< 1 µg/10 <sup>6</sup> cells
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG20702 Rat IgG2a Isotype Control antibody [KLH/G2a-1-1] (FITC)  
FLISA image

FLISA: The plate was coated with [ARG21694](#) Mouse anti-Rat IgG1 antibody [G17E7], [ARG21699](#) Mouse anti-Rat IgG2a antibody [2A8F4], [ARG21705](#) Mouse anti-Rat IgG2b antibody [2B10A8], [ARG21711](#) Mouse anti-Rat IgG2c antibody [2C8F1] and [ARG23828](#) Mouse anti-Rat IgM antibody [M2A1]. Serially diluted ARG20702 Rat IgG2a Isotype Control antibody [KLH/G2a-1-1] (FITC) was captured and fluorescence intensity quantified.