

## ARG20709 Rat IgG2a Isotype Control antibody [KLH/G2a-1-1] (PE-Cyanine 7)

Package: 50 µg  
Store at: 4°C

### Summary

Product Description	PE-Cyanine 7-conjugated Rat Monoclonal antibody [KLH/G2a-1-1] as a negative control antibody for Rat IgG2a
Tested Application	BL, ELISA, FACS, FLISA, IHC-Fr, IHC-P
Specificity	KLH
Host	Rat
Clonality	Monoclonal
Clone	KLH/G2a-1-1
Isotype	IgG2a, kappa
Target Name	Rat IgG2a
Conjugation	PE-Cyanine 7

### Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	ELISA	Assay-dependent
	FACS	< 0.2 µg/10 <sup>6</sup> cells
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.