

## ARG20724

### Rat IgG2b Isotype Control antibody [KLH/G2b-1-2] (APC-Cyanine 7)

Package: 50 µg

Store at: 4°C

## Summary

Product Description	APC-Cyanine 7-conjugated Rat Monoclonal antibody [KLH/G2b-1-2] as a negative control antibody for Rat IgG2b
Tested Application	ELISA, FACS, FLISA, IHC-Fr, IHC-P
Specificity	KLH
Host	Rat
Clonality	Monoclonal
Clone	KLH/G2b-1-2
Isotype	IgG2b, kappa
Target Name	Rat IgG2b
Conjugation	APC-Cyanine 7

## Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	< 0.1 µg/10 <sup>6</sup> cells
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

## Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.