

ARG20727 Rat IgG2c Isotype Control antibody [SB68b] (Biotin)

Package: 100 µg
Store at: 4°C

Summary


Product Description	Biotin-conjugated Rat Monoclonal antibody [SB68b] as a negative control antibody for Rat IgG2c
Tested Application	ELISA, FACS, FLISA
Specificity	KLH
Host	Rat
Clonality	Monoclonal
Clone	SB68b
Isotype	IgG2c, kappa
Target Name	Rat IgG2c
Conjugation	Biotin

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	< 1 µg/10 ⁶ cells
	FLISA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG20727 Rat IgG2c Isotype Control antibody [SB68b] (Biotin) ELISA image

ELISA: The plate was coated with Mouse anti-Rat IgG1 antibody, [ARG21699](#) Mouse anti-Rat IgG2a antibody [2A8F4], [ARG21705](#) Mouse anti-Rat IgG2b antibody [2B10A8], [ARG21711](#) Mouse anti-Rat IgG2c antibody [2C8F1], and [ARG23828](#) Mouse anti-Rat IgM antibody [M2A1]. Serially diluted ARG20727 Rat IgG2c Isotype Control antibody [SB68b] (Biotin) was captured followed by [ARG23912](#) Streptavidin (HRP) and quantified.