

## ARG20729 Rat IgG2c Isotype Control antibody [SB68b] (APC)

Package: 50 µg  
Store at: 4°C

### Summary


Product Description	APC-conjugated Rat Monoclonal antibody [SB68b] as a negative control antibody for Rat IgG2c
Tested Application	FACS, FLISA
Specificity	KLH
Host	Rat
Clonality	Monoclonal
Clone	SB68b
Isotype	IgG2c, kappa
Target Name	Rat IgG2c
Conjugation	APC

### Application Instructions

Application table	Application	Dilution
	FACS	< 0.1 µg/10 <sup>6</sup> cells
	FLISA	< 1 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG20729 Rat IgG2c Isotype Control antibody [SB68b] (APC) FLISA image

FLISA: The plate was coated with Mouse anti-Rat IgG1 antibody, [ARG21699](#) Mouse anti-Rat IgG2a antibody [2A8F4], [ARG21705](#) Mouse anti-Rat IgG2b antibody [2B10A8], [ARG21711](#) Mouse anti-Rat IgG2c antibody [2C8F1] and [ARG23828](#) Mouse anti-Rat IgM antibody [M2A1]. Serially diluted ARG20729 Rat IgG2c Isotype Control antibody [SB68b] (APC) was captured and fluorescence intensity quantified.