

ARG20739 Syrian Hamster IgG Isotype Control antibody

Package: 100 µg

Store at: -20°C

Summary

Product Description	Syrian Hamster Polyclonal antibody as a negative control antibody for Hamster IgG
Tested Application	ELISA, FACS, FLISA
Host	Syrian Hamster
Clonality	Polyclonal
Isotype	IgG
Target Name	Syrian Hamster IgG
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	Assay-dependent
	FLISA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.