

**ARG20743**  
**Syrian Hamster IgG Isotype Control antibody (APC)**Package: 50 µg  
Store at: 4°C

### Summary

Product Description	APC-conjugated Syrian Hamster Polyclonal antibody as a negative control antibody for Hamster IgG
Tested Application	ELISA, FACS, FLISA
Host	Syrian Hamster
Clonality	Polyclonal
Isotype	IgG
Target Name	Syrian Hamster IgG
Conjugation	APC

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	< 0.2 µg/10 <sup>6</sup> cells
	FLISA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.