

ARG20748 Syrian Hamster IgG Isotype Control antibody (APC-Cyanine 5.5)

Package: 50 µg
Store at: 4°C

Summary

Product Description	APC-Cyanine 5.5-conjugated Syrian Hamster Polyclonal antibody as a negative control antibody for Hamster IgG
Tested Application	ELISA, FACS, FLISA
Host	Syrian Hamster
Clonality	Polyclonal
Isotype	IgG
Target Name	Syrian Hamster IgG
Conjugation	APC-Cyanine 5.5

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	< 0.2 µg/10 ⁶ cells
	FLISA	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.