

ARG20753 Goat IgG Fc Isotype Control antibody

Package: 1 mg
Store at: -20°C

Summary


Product Description	Goat Polyclonal antibody as a negative control antibody for Goat IgG Fc
Tested Application	ELISA
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	Goat IgG Fc
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	< 2 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified
Buffer	BBS (pH 8.2)
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG20753 Goat IgG Fc Isotype Control antibody ELISA image

ELISA: The plate was coated with serially diluted ARG20753 Goat IgG Fc Isotype Control antibody. Immunoglobulin was detected with [ARG23853](#) Mouse anti-Goat IgG Fc antibody [SB115h] (Biotin) and Rabbit anti-Goat IgG (Fab) antibody (Biotin) followed by [ARG23912](#) Streptavidin (HRP) and quantified.