

## ARG20755 Sheep IgG Isotype Control antibody

Package: 1 mg  
Store at: -20°C

### Summary

Product Description	Sheep Polyclonal antibody as a negative control antibody for Sheep IgG
Tested Application	ELISA, FACS, FLISA
Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Target Name	Sheep IgG
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	< 1 µg/10 <sup>6</sup> cells
	FLISA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Purified
Buffer	BBS (pH 8.2)
Concentration	5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.