

## ARG20758 Rhesus Monkey IgG Isotype Control antibody

Package: 500 µg  
Store at: -20°C

### Summary

Product Description	Rhesus Monkey Polyclonal antibody as a negative control antibody for Rhesus Monkey IgG
Tested Application	ELISA
Host	Rhesus Monkey
Clonality	Polyclonal
Isotype	IgG
Target Name	Rhesus Monkey IgG
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	< 2 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Purified
Buffer	BBS (pH 8.2)
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.