

ARG20762 Human IgG Isotype Control antibody (HRP)

Package: 500 µl
Store at: 4°C

Summary


Product Description	HRP-conjugated Human Polyclonal antibody as a negative control antibody for Human IgG
Tested Application	ELISA
Host	Human
Clonality	Polyclonal
Isotype	IgG
Target Name	Human IgG
Conjugation	HRP

Application Instructions

Application table	Application	Dilution
	ELISA	1:2000 - 1:8000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	50% PBS (pH 7.4) and 50% Glycerol
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG20762 Human IgG Isotype Control antibody (HRP) ELISA image

ELISA: The plate was coated with [ARG23878](#) Goat anti-Human IgG antibody (pre-adsorbed), [ARG21644](#) Goat anti-Human IgM antibody (pre-adsorbed), and [ARG21909](#) Goat anti-Human IgA antibody (pre-adsorbed). Serially diluted ARG20762 Human IgG Isotype Control antibody (HRP) was captured and quantified.