

ARG20767 Human IgG1 Kappa Isotype Control antibody (low endotoxin)

Package: 100 µg
Store at: -20°C

Summary

Product Description	Azide free and low endotoxin Human Polyclonal antibody as a negative control antibody for Human IgG1 Kappa
Tested Application	ELISA, FACS, FLISA, IP
Host	Human
Clonality	Polyclonal
Isotype	IgG1, kappa
Target Name	Human IgG1 Kappa
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	< 2 µg/ml
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IP	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Purified
Purification Note	Low endotoxin
Buffer	PBS
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.