

ARG20773 Human IgA Kappa Isotype Control antibody (AP)

Package: 500 µl
Store at: 4°C

Summary

Product Description	AP-conjugated Human Polyclonal antibody as a negative control antibody for Human IgA Kappa
Tested Application	ELISA
Host	Human
Clonality	Polyclonal
Isotype	IgA, kappa
Target Name	Human IgA Kappa
Conjugation	AP

Application Instructions

Application table	Application	Dilution
	ELISA	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	50 mM Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl ₂ and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.