

**ARG20776**  
**Human IgA Lambda Isotype Control antibody (AP)**Package: 500 µl  
Store at: 4°C

### Summary

Product Description	AP-conjugated Human Polyclonal antibody as a negative control antibody for Human IgA Lambda
Tested Application	ELISA
Host	Human
Clonality	Polyclonal
Isotype	IgA, lambda
Target Name	Human IgA Lambda
Conjugation	AP

### Application Instructions

Application table	Application	Dilution
	ELISA	1:1000 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	50 mM Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl <sub>2</sub> and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.