

ARG21085 anti-MHC Class I H2 Dk antibody [15-5-5S] (FITC)

Package: 100 µg
Store at: 4°C

Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [15-5-5S] recognizes MHC Class I H2 Dk
Tested Reactivity	Ms
Tested Application	BL, ELISA, FACS
Specificity	Mouse H-2Dk. The clone 15-5-5 reacts with the H-2Dk class I alloantigen. It cross-reacts with H-2Kd and with cells from mice with the H-2f haplotype. The antibody does not react with other (e.g., a, b, p, q, r, s) haplotypes.
Host	Mouse
Clonality	Monoclonal
Clone	15-5-5S
Isotype	IgG2a, kappa
Target Name	MHC Class I H2 Dk
Species	Mouse
Immunogen	C3H mouse splenocytes
Conjugation	FITC
Alternate Names	H-2 class I histocompatibility antigen, D-B alpha chain; H2-D; B; H-2D

Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	ELISA	Assay-dependent
	FACS	< 1 µg/10 ⁶ cells
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Database links	GeneID: 14964 Mouse Swiss-port # P01899 Mouse
Gene Symbol	H2-D1
Gene Full Name	histocompatibility 2, D region locus 1
Function	Involved in the presentation of foreign antigens to the immune system. [UniProt]
Calculated Mw	41 kDa