


## ARG21132 anti-CD172a / SIRP alpha antibody [74-22-15] (Biotin)

Package: 250 µg  
Store at: 4°C

### Summary

Product Description	Biotin-conjugated Mouse Monoclonal antibody [74-22-15] recognizes CD172a / SIRP alpha
Tested Reactivity	Pig
Tested Application	ELISA, FACS, ICC/IF, IHC-Fr, IHC-P
Specificity	

### Images



ARG21132 anti-CD172a / SIRP alpha antibody [74-22-15] (Biotin)  
FACS image

Flow Cytometry: Porcine small lymphocytes, monocytes, and granulocytes stained with [ARG21132](#) anti-CD172a / SIRP alpha antibody [74-22-15] (Biotin) and [ARG21114](#) anti-CD3e antibody [PPT3] (PE).