


ARG21133 anti-CD172a / SIRP alpha antibody [74-22-15] (PE)

Package: 50 µg
Store at: 4°C

Summary

Product Description	PE-conjugated Mouse Monoclonal antibody [74-22-15] recognizes CD172a / SIRP alpha
Tested Reactivity	Pig
Tested Application	FACS, ICC/IF, IHC-Fr, IHC-P
Specificity	

Images



ARG21133 anti-CD172a / SIRP alpha antibody [74-22-15] (PE) FACS image

Flow Cytometry: Porcine peripheral blood granulocytes stained with ARG21133 anti-CD172a / SIRP alpha antibody [74-22-15] (PE).