

ARG21203 anti-MHC Class I antibody [F21-2] (PE)

Package: 50 µg
Store at: 4°C

Summary

Product Description	PE-conjugated Mouse Monoclonal antibody [F21-2] recognizes MHC Class I
Tested Reactivity	Chk, Turkey
Tested Application	BL, EM, FACS, IHC-Fr, Puri, WB
Specificity	Chicken/Turkey MHC Class I. The clone F21-2 also reacts with turkey MHC Class I.
Host	Mouse
Clonality	Monoclonal
Clone	F21-2
Isotype	IgG1, kappa
Target Name	MHC Class I
Antigen Species	Chicken
Conjugation	PE

Application Instructions

Application table	Application	Dilution
	BL	Assay-dependent
	EM	Assay-dependent
	FACS	< 0.2 µg/10 ⁶ cells
	IHC-Fr	Assay-dependent
	Puri	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.