

## ARG21212 anti-Bu1a antibody [21-1A4]

Package: 250 µg  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [21-1A4] recognizes Bu1a
Tested Reactivity	Chk
Tested Application	FACS, ICC/IF, IHC-Fr, IP
Specificity	Chicken Bu-1a. The clone 21-1A4 does not react with cells from CHA and H.B14ab strains by immunofluorescence.
Host	Mouse
Clonality	Monoclonal
Clone	21-1A4
Isotype	IgG1, kappa
Target Name	Bu1a
Species	Chicken
Immunogen	Bursa cells from one-day-old H.B15 (Bu-1a/b) strain chickens
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	FACS	< 1 µg/10 <sup>6</sup> cells
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent
	IP	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.