

ARG21213 anti-Bu1a antibody [21-1A4] (FITC)

Package: 100 µg
Store at: 4°C

Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [21-1A4] recognizes Bu1a
Tested Reactivity	Chk
Tested Application	FACS, ICC/IF, IHC-Fr
Specificity	Chicken Bu-1a. The clone 21-1A4 does not react with cells from CHA and H.B14ab strains by immunofluorescence.
Host	Mouse
Clonality	Monoclonal
Clone	21-1A4
Isotype	IgG1, kappa
Target Name	Bu1a
Antigen Species	Chicken
Immunogen	Bursa cells from one-day-old H.B15 (Bu-1a/b) strain chickens
Conjugation	FITC

Application Instructions

Application table	Application	Dilution
	FACS	< 1 µg/10 ⁶ cells
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.