

ARG21596
Rat anti-Mouse IgE antibody [23G3] (low endotoxin)Package: 250 µg
Store at: -20°C

Summary

Product Description	Azide free and low endotoxin Rat Monoclonal antibody [23G3] recognizes Mouse IgE
Tested Reactivity	Ms
Tested Application	ELISA, ELISPOT, FACS, FLISA, IHC-Fr, IHC-Wmt, Neut, WB
Host	Rat
Clonality	Monoclonal
Clone	23G3
Isotype	IgG1, kappa
Target Name	IgE
Species	Mouse
Immunogen	Monoclonal anti-DNP mouse IgE
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	ELISPOT	Assay-dependent
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-Wmt	Assay-dependent
	Neut	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification Note	Low endotoxin
Buffer	PBS
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot

and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.