

ARG21679

Goat anti-Human Kappa Light Chain antibody, F(ab')₂ fragment (HRP),
pre-adsorbed

Package: 500 µl

Store at: 4°C

Summary

Product Description	HRP-conjugated Goat Polyclonal antibody recognizes Human Kappa Light Chain
Tested Reactivity	Hu
Tested Application	ELISA, ELISPOT, FACS, FLISA, ICC/IF, WB
Specificity	The antibody reacts with Human kappa light chains. The antibody is pre-adsorbed with Human lambda light chains, so the antibody may not react with Human lambda light chains, but may react with kappa light chains from other species.
Host	Goat
Clonality	Polyclonal
Isotype	F(ab') ₂ IgG
Target Name	Kappa Light Chain
Antigen Species	Human
Conjugation	HRP

Application Instructions

Pre Adsorbed	Human lambda light chains.														
Application table	<table><thead><tr><th>Application</th><th>Dilution</th></tr></thead><tbody><tr><td>ELISA</td><td>1:4000 - 1:8000</td></tr><tr><td>ELISPOT</td><td>Assay-dependent</td></tr><tr><td>FACS</td><td>Assay-dependent</td></tr><tr><td>FLISA</td><td>Assay-dependent</td></tr><tr><td>ICC/IF</td><td>Assay-dependent</td></tr><tr><td>WB</td><td>Assay-dependent</td></tr></tbody></table>	Application	Dilution	ELISA	1:4000 - 1:8000	ELISPOT	Assay-dependent	FACS	Assay-dependent	FLISA	Assay-dependent	ICC/IF	Assay-dependent	WB	Assay-dependent
Application	Dilution														
ELISA	1:4000 - 1:8000														
ELISPOT	Assay-dependent														
FACS	Assay-dependent														
FLISA	Assay-dependent														
ICC/IF	Assay-dependent														
WB	Assay-dependent														
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.														

Properties


Form	Liquid
Buffer	50% PBS (pH 7.4) and 50% Glycerol
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be

gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.

Images



ARG21679 Goat anti-Human Kappa Light Chain antibody, F(ab')₂ fragment (HRP) (pre-adsorbed) ELISA image

ELISA: The plate was coated with purified Human IgGκ, IgMκ, IgAκ, IgGλ, IgMλ, and IgAλ. Immunoglobulins were detected with serially diluted ARG21679 Goat anti-Human Kappa Light Chain antibody, F(ab')₂ fragment (HRP) (pre-adsorbed).