

**ARG21747**  
**Goat anti-Chicken IgY (H+L) antibody (HRP)**Package: 500 µl  
Store at: 4°C

### Summary

Product Description	HRP-conjugated Goat Polyclonal antibody recognizes Chicken IgY (H+L)
Tested Reactivity	Chk
Tested Application	ELISA, FACS, FLISA, ICC/IF, IHC-P, WB
Specificity	Reacts with the heavy and light chains of chicken IgY
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgY (H+L)
Antigen Species	Chicken
Conjugation	HRP

### Application Instructions

Application table	Application	Dilution
	ELISA	1:4000 - 1:8000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	50% PBS (pH 7.4) and 50% Glycerol
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.