

## ARG21793 Mouse anti-Chicken IgA antibody [A-1]

Package: 250 µg  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [A-1] recognizes Chicken IgA
Tested Reactivity	Chk
Tested Application	ELISA, FACS, ICC/IF, IHC-Fr, IHC-P, IP
Host	Mouse
Clonality	Monoclonal
Clone	A-1
Isotype	IgG2b, kappa
Target Name	IgA
Species	Chicken
Immunogen	Affinity purified chicken Ig or isolated lymphocytes
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	Assay-dependent
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	IP	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.