

## ARG21796 Mouse anti-Chicken IgA antibody [A-1] (PE)

Package: 100 µg  
Store at: 4°C

### Summary

Product Description	PE-conjugated Mouse Monoclonal antibody [A-1] recognizes Chicken IgA
Tested Reactivity	Chk
Tested Application	ELISA, FACS, ICC/IF, IHC-Fr, IHC-P
Host	Mouse
Clonality	Monoclonal
Clone	A-1
Isotype	IgG2b, kappa
Target Name	IgA
Species	Chicken
Immunogen	Affinity purified chicken Ig or isolated lymphocytes
Conjugation	PE

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	< 0.2 µg/10 <sup>6</sup> cells
	ICC/IF	< 10 µg/ml
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

**Note**

For laboratory research only, not for drug, diagnostic or other use.