

ARG21797 Mouse anti-Chicken Lambda Light Chain antibody [L-1]

Package: 250 µg

Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [L-1] recognizes Chicken Lambda Light Chain
Tested Reactivity	Chk
Tested Application	FACS, ICC/IF, IHC-Fr, IP, WB
Host	Mouse
Clonality	Monoclonal
Clone	L-1
Isotype	IgG1, kappa
Target Name	Lambda Light Chain
Species	Chicken
Immunogen	Affinity purified chicken Ig or isolated lymphocytes
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	FACS	Assay-dependent
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent
	IP	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.