

ARG21820 Mouse anti-Human IgG1 (Hinge) antibody [4E3] (FITC)

Package: 250 µg
Store at: 4°C

Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [4E3] recognizes Human IgG1 (Hinge)
Tested Reactivity	Hu
Tested Application	FACS, FLISA, IHC-P, WB
Host	Mouse
Clonality	Monoclonal
Clone	4E3
Isotype	IgG1, kappa
Target Name	IgG1 (Hinge)
Antigen Species	Human
Conjugation	FITC


Application Instructions

Application table	Application	Dilution
	FACS	Assay-dependent
	FLISA	1:200 - 1:400
	IHC-P	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG21820 Mouse anti-Human IgG1 (Hinge) antibody [4E3] (FITC)
FLISA image

FLISA: The plate was coated with purified Human IgG1, IgG2, IgG3 and IgG4. Immunoglobulins were detected with serially diluted ARG21820 Mouse anti-Human IgG1 (Hinge) antibody [4E3] (FITC).