

ARG21831 Mouse anti-Human IgG2 (Fc) antibody [HP6002] (PE)

Package: 100 µg
Store at: 4°C

Summary

Product Description	PE-conjugated Mouse Monoclonal antibody [HP6002] recognizes Human IgG2 (Fc)
Tested Reactivity	Hu
Tested Application	FACS, FLISA, IHC-Fr, WB
Host	Mouse
Clonality	Monoclonal
Clone	HP6002
Isotype	IgG1, kappa
Target Name	IgG2 (Fc)
Species	Human
Immunogen	Human IgG2 myeloma protein
Conjugation	PE


Application Instructions

Application table	Application	Dilution
	FACS	Assay-dependent
	FLISA	< 1 ug/ml
	IHC-Fr	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG21831 Mouse anti-Human IgG2 (Fc) antibody [HP6002] (PE)
FLISA image

FLISA: The plate was coated with purified Human IgG1, IgG2, IgG3 and IgG4. Immunoglobulins were detected with serially diluted ARG21831 Mouse anti-Human IgG2 (Fc) antibody [HP6002] (PE).